

```

1  /*
2  Object-Location Task
3
4  %%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%
5  %%%%%%%%% Version History %%%%%%%%%
6  %%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%
7  - 06-10-2019; created by Chris Gill: created the backbone of the experiment
8  - 10-19-2019; edited by Alex He: changed the structure of the task, improved the
  iteration code block
9  - 11-09-2020; edited by Alex He: added triggers, added version history, added
  22/23 start/end triggers, updated header parameters, enabled proper logfile saving
10 - 11-11-2020; edited by Alex He: changed pulse width from 10ms to 5ms, added more
  wait_delays since response_port_output is true
11 - 01-19-2021; edited by Alex He: updated logfile naming
12 - 10-05-2022; edited by Anthony Edgar: changed task to be compatible with c-pod.
13 - 03-05-2024; edited by Alex He: added fixation cross to the 8-s blank screen
14 %%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%
15
16 *Note: After each round of practice trials the experimenter will be asked to
  "Please press keys to continue."
17 In order to continue to the trial press 1, in order to repeat practice trials
  press 2.
18
19 */
20 #####
  #####
21 # Header
22 scenario = "Object_Location_Task";
23 active_buttons = 5;
24 button_codes = 11,12,13,14,15; # keys: 1 = spacebar, 2 = 'Y' (Yes), 3 = 'N' (No),
  4 = #1 (continue), 5 = #2 (repeat practice)
25 default_font = "Helvetica";
26 default_background_color = 255,255,255;
27 default_text_color = 0,0,0;
28 write_codes = true;
29 response_port_output = true;
30 default_output_port = 1;
31 pulse_width = 5;
32 #####
  #####
33 #SDL
34 begin;
35
36 #load image array
37 array {
38 LOOP $i 179;
39 $k = '$i+1';
40 bitmap { filename = "$k.jpg"; description = "Image$k"; width = 320; height = 225;
  }"image$k";
41 ENDLLOOP;
42 } pictures;
43
44 array {
45 LOOP $i 79;
46 $k = '$i+180';
47 bitmap { filename = "$k.jpg"; description = "Image$k"; width = 225; height = 320;
  }"image$k";
48 ENDLLOOP;
49 } pictures_b;
50
51 #circle object
52 bitmap { filename = "circle_obj.jpg"; description = "Circle"; width = 200; height
  = 200; }"circle";
53

```

```

54 #Grid
55 line_graphic { # x1, y1, x2, y2
56     coordinates = -500, -500, -500, 500; #outer vert left
57     coordinates = 500, -500, 500, 500; #outer vert right
58     coordinates = -500, 498, 500, 498; #outer hor top
59     coordinates = -500, -497, 500, -497; #outer hor bottom
60     coordinates = -166.66, -500, -166.66, 500; #inner vert left
61     coordinates = 166.66, -500, 166.66, 500; #inner vert right
62     coordinates = -500, 166.66, 500, 166.66; #inner hor top
63     coordinates = -500, -166.66, 500, -166.66; #inner hor bottom
64     line_width = 6;
65     line_color = 0, 0, 0, 255;
66 }grid;
67
68 #Prompt screen
69 text {font_size = 60; font_color = 0,0,0; background_color = 255,255,255; caption
= "Remember Object";} Obj_text;
70 text {font_size = 60; font_color = 0,0,0; background_color = 255,255,255; caption
= "Remember Location";} Loc_text;
71 text {font_size = 60; font_color = 0,0,0; background_color = 255,255,255; caption
= "Remember Object + Location";} Both_text;
72
73 #Blank screen
74 box {height = 1920; width = 1080; color = 255, 255,255; } blank_box;
75 text { caption = "+"; font_size = 96; font_color = 0,0,0; } Crosshair_text;
76
77 #####
78 #Instructions + practice pics
79 picture {
80     text {
81         caption = "<u>INSTRUCTIONS:</u>";
82         font_size = 44;
83         formatted_text = true;
84     };
85     x = 0;y = 380;
86
87     text { caption =
88         "In this experiment you will see a grid with 9 squares.
89         3 objects will appear one at a time at different locations
90         within the grid. You will be asked to look at these objects,
91         then after a short delay you will be tested on how well you
92         can remember them.";
93         font_size = 44;
94         text_align = align_left;
95     } Instruction1;
96     x = 0;y = 10;
97
98     text {
99         caption = "Press any key to continue.";
100         font_size = 36;
101     };
102     x = 0;y = -500;
103 } instructions1;
104
105
106 picture {
107     text {
108         caption = "<u>INSTRUCTIONS:</u>";
109         font_size = 44;
110         formatted_text = true;
111     };
112     x = 0;y = 380;
113

```

```

114     text { caption =
115         "There will be two different types of trials for this part of the
116         experiment. The first type is called a <b>'Remember Object'</b>
117         trial. The second is called a <b>'Remember Location'</b> trial. We
118         will now explain each type of trial separately.
119
120         While we explain, you will be shown some pictures of the trials
121         and afterwards you will have a chance to do some practice trials.";
122         font_size = 44;
123         text_align = align_left;
124         formatted_text = true;
125     } Instruction2;
126     x = 0;y = 10;
127
128     text {
129         caption = "Press any key to continue.";
130         font_size = 36;
131     };
132     x = 0;y = -500;
133 } instructions2;
134
135
136 picture {
137     text { caption =
138         "<u><b>Remember Object trial:</b></u> In this type of trial you need to
remember
139         the identity of the objects shown to you.
140
141         You will see 3 objects appearing one at a time followed by an
142         8 second delay. You will then be shown a prompt screen saying
143         <b>Remember Object</b>. This prompt screen will be followed by a test
144         object <u>in the center of the grid</u>.
145
146         'Remember Object' tells you that you need to decide whether this
147         test object was one of the 3-object sequence just shown to you. ";
148         font_size = 44;
149         text_align = align_left;
150         formatted_text = true;
151     } Instruction3;
152     x = 0;y = 0;
153
154     text {
155         caption = "Press any key to continue.";
156         font_size = 32;
157     };
158     x = 0;y = -500;
159 } instructions3;
160
161
162 picture {
163     text { caption =
164         "<u><b>Remember Object trial:</b></u>";
165         font_size = 44;
166         formatted_text = true;
167     };
168     x = 0;y = 450;
169
170     bitmap { filename = "Obj_trial_instructions_pic.jpg"; width = 1650; height =
751; };
171     x = 0;y = -20;
172
173     text {
174         caption = "Press any key to continue.";
175         font_size = 32;

```

```

176     };
177     x = 0;y = -500;
178 } instructions4;
179
180
181 picture {
182     bitmap { filename = "Obj_trial_instructions_pic.jpg"; width = 1581; height =
183 720; };
184     x = 0;y = -160;
185     text { caption =
186         "<u><b>Remember Object trial:</b></u>"
187         "If you recognize the test object, press the <b>'Y'</b> key to indicate
188         <b>YES</b> the test object was one of
189         the 3-object sequence.
190         If you do not recognize the test object, press <b>'N'</b> key to indicate
191         <b>NO</b> the test object was
192         <u>NOT</u> one of the 3-object sequence.";
193         font_size = 34;
194         text_align = align_left;
195         formatted_text = true;
196     } Instruction5;
197     x = 0;y = 360;
198     text {
199         caption = "Press any key to continue.";
200         font_size = 32;
201     };
202     x = 0;y = -500;
203 } instructions5;
204
205
206 picture {
207     text { caption =
208         "<u><b>Remember Location trial:</b></u> In this type of trial you need to
209 remember
210         the location of the objects shown to you.
211         You will see 3 objects appearing one at a time followed by an
212         8 second delay. You will then be shown a prompt screen saying
213         <b>Remember Location</b>. This prompt screen will be followed by a
214         dot in one of the squares of the grid.
215         'Remember Location' tells you that you need to decide whether
216         this square was previously occupied by <u>any object</u> in the
217         3-object sequence.";
218         font_size = 44;
219         text_align = align_left;
220         formatted_text = true;
221     } Instruction6;
222     x = 0;y = 0;
223     text {
224         caption = "Press any key to continue.";
225         font_size = 32;
226     };
227     x = 0;y = -500;
228 } instructions6;
229
230
231
232
233 picture {
234     text { caption =
235         "<u><b>Remember Location trial:</b></u>";

```

```

236     font_size = 44;
237     formatted_text = true;
238 };
239 x = 0;y = 450;
240
241 bitmap { filename = "Loc_trial_instructions_pic.jpg"; width = 1650; height =
732; };
242 x = 0;y = -20;
243
244 text {
245     caption = "Press any key to continue.";
246     font_size = 32;
247 };
248 x = 0;y = -500;
249 } instructions7;
250
251
252 picture {
253     bitmap { filename = "Loc_trial_instructions_pic.jpg"; width = 1585; height =
704; };
254     x = 0;y = -160;
255
256     text { caption =
257         "<u><b>Remember Location trial:</b></u>"
258         "If the square indicated by the dot was previously occupied, press the"
259         "<b>'Y'</b> key to indicate <b>YES</b>"
260         "one of the objects in the 3-object sequence was in this square."
261         "If the square indicated by the dot was <u>NOT</u> previously occupied,"
262         "press <b>'N'</b> key to indicate <b>NO</b>"
263         "none of the objects in the 3-object sequence were in this square.";
264     font_size = 34;
265     text_align = align_left;
266     formatted_text = true;
267 } Instruction8;
268 x = 0;y = 360;
269
270 text {
271     caption = "Press any key to continue.";
272     font_size = 32;
273 };
274 x = 0;y = -500;
275 } instructions8;
276
277 picture {
278     text {
279         caption = "<u>INSTRUCTIONS:</u>";
280         font_size = 44;
281         formatted_text = true;
282     };
283     x = 0;y = 380;
284
285     text { caption =
286         "<u><b>Review:</b></u>"
287         "There are 2 types of trials in this part of the experiment:"
288         "<u><b>Remember Object</b></u> and <u><b>Remember Location</b></u>."
289
290         "Note that Remember Object and Remember Location trials will"
291         "be intermixed throughout this part of the experiment. This"
292         "means that you will not know whether you need to respond to"
293         "the Object or the Location until you see the Prompt screen." ;
294         font_size = 44;

```

```

296     text_align = align_left;
297     formatted_text = true;
298 } Instruction9;
299 x = 0;y = 10;
300
301 text {
302     caption = "Click the space bar to start the practice trials.";
303     font_size = 36;
304 };
305 x = 0;y = -500;
306 } instructions9;
307
308
309 picture {
310     text {
311         caption = "<u>INSTRUCTIONS:</u>";
312         font_size = 44;
313         formatted_text = true;
314     };
315     x = 0;y = 380;
316
317     text { caption =
318         "In this second part of the experiment you will again see a
319         grid with 9 squares. 3 objects will appear one at a time at
320         different locations within the grid. You will be asked to look
321         at these objects, and after a short delay, you will be tested
322         on how well you can remember <b>both Objects and their Locations</b>.";
323         font_size = 44;
324         text_align = align_left;
325         formatted_text = true;
326     } Instruction10;
327     x = 0;y = 10;
328
329     text {
330         caption = "Press any key to continue.";
331         font_size = 36;
332     };
333     x = 0;y = -500;
334 } instructions10;
335
336
337 picture {
338     text { caption =
339         "<u><b>Remember Object and Location trial:</b></u> In this type of trial
340         you need
341         to remember both the identity of the objects and their locations.
342
343         As before, you will see 3 objects appearing one at a time followed by
344         an 8 second delay. You will then be shown a prompt screen saying
345         <b>Remember Object and Location</b>. This prompt screen will be
346         followed by a test object <u>in one of the squares of the grid</u>.
347
348         'Remember Object and Location' tells you that you need to decide
349         whether this object is in the right location. What this means is that
350         the test object is in the same square as it was in the 3-object sequence.";
351         font_size = 44;
352         text_align = align_left;
353         formatted_text = true;
354     } Instruction11;
355     x = 0;y = 0;
356
357     text {
358         caption = "Press any key to continue.";
359         font_size = 32;

```

```

359     };
360     x = 0;y = -500;
361 } instructions11;
362
363
364 picture {
365     text { caption =
366         "<u><b>Remember Object and Location trial:</b></u>";
367         font_size = 44;
368         formatted_text = true;
369     };
370     x = 0;y = 450;
371
372     bitmap { filename = "ObjLoc_trial_instructions_pic.jpg"; width = 1650; height =
373         727; };
374     x = 0;y = -20;
375
376     text {
377         caption = "Press any key to continue.";
378         font_size = 32;
379     };
380     x = 0;y = -500;
381 } instructions12;
382
383 picture {
384     bitmap { filename = "ObjLoc_trial_instructions_pic.jpg"; width = 1358; height =
385         599; };
386     x = 0;y = -220;
387
388     text { caption =
389         "<u><b>Remember Object and Location trial:</b></u>
390         If the test object is in the right location, press the <b>'Y'</b> key to
391         indicate <b>YES</b> that the test object
392         is in the same square as it was in the 3-object sequence.
393
394         If the test object is <u>NOT</u> in the right location, press the
395         <b>'N'</b> key to indicate <b>NO</b> the test object
396         is <u>NOT</u> in the same square as it was in the 3-object sequence.";
397         font_size = 34;
398         text_align = align_left;
399         formatted_text = true;
400     } Instruction13;
401     x = 0;y = 360;
402
403     text {
404         caption = "Click the space bar to start the practice trials.";
405         font_size = 32;
406     };
407     x = 0;y = -500;
408 } instructions13;
409
410 #####
411 #####
412 #Trials List
413
414 #Instruction
415 trial {
416     trial_type = first_response;
417     trial_duration = forever;
418     stimulus_event {
419         picture instructions1;
420         code = "Instructions";
421     } Instructions_event;

```

```

418 }Instructions_trial;
419
420 #Blank grid
421 trial {
422     trial_type = fixed;
423     trial_duration = 500;
424     picture {
425         line_graphic grid;
426         x = 0; y = 0;
427     };
428 } blank_grid_trial;
429
430 #Blank screen for 8s delay
431 trial {
432     trial_type = fixed;
433     trial_duration = 8000;#####Change back to 8000 after testing
434     picture {
435         box blank_box;
436         x = 0; y = 0;
437
438         text Crosshair_text;
439         x = 0; y = 0;
440     };
441     port_code = 2;
442     code = "8s Blank";
443 } wait_delay;
444
445 #Remember Object prompt screen
446 trial {
447     trial_type = fixed;
448     trial_duration = 1000;
449     picture {
450         background_color = 255,255,255;
451         text Obj_text;
452         x = 0; y = 0;
453     };
454     port_code = 25;
455     code = "Remember Object Prompt";
456 } Obj_Condition_Message;
457
458 #Remember Location prompt screen
459 trial {
460     trial_type = fixed;
461     trial_duration = 1000;
462     picture {
463         background_color = 255,255,255;
464         text Loc_text;
465         x = 0; y = 0;
466     };
467     port_code = 26;
468     code = "Remember Location Prompt";
469 } Loc_Condition_Message;
470
471 #Remember Object and Location prompt screen
472 trial {
473     trial_type = fixed;
474     trial_duration = 1000;
475     picture {
476         background_color = 255,255,255;
477         text Both_text;
478         x = 0; y = 0;
479     };
480     port_code = 27;
481     code = "Remember Obj+Loc Prompt";

```



```

482 } Conjunction_Condition_Message;
483
484 #Practice Encoding trial
485 trial {
486     trial_type = fixed;
487     trial_duration = 1000;
488     stimulus_event {
489         picture {
490             background_color = 255,255,255;
491             line_graphic grid;
492             x = 0; y = 0;
493
494             bitmap image1;
495             x = 0; y = 0;
496         } practice_pic;
497         port_code = 7;
498         code = "Practice Encoding";
499     };
500 } practice_study_trial;
501
502 #Encoding trial
503 trial {
504     trial_type = fixed;
505     trial_duration = 1000;
506     stimulus_event {
507         picture {
508             background_color = 255,255,255;
509             line_graphic grid;
510             x = 0; y = 0;
511
512             bitmap image1;
513             x = 0; y = 0;
514         } pic;
515         port_code = 3;
516         code = "Encoding Trial";
517     };
518 } study_trial;
519
520 #Practice Remember Object test trial
521 trial {
522     trial_type = correct_response;
523     trial_duration = forever;
524     stimulus_event {
525         picture {
526             background_color = 255,255,255;
527             line_graphic grid;
528             x = 0; y = 0;
529
530             bitmap image1;
531             x = 0; y = 0;
532         } practice_Obj_test_pic;
533         target_button = 2,3;
534         port_code = 8;
535         code = "Practice Remember Object Test";
536     };
537 } practice_Obj_test_trial;
538
539 #Remember Object test trial
540 trial {
541     trial_type = correct_response;
542     trial_duration = forever;
543     stimulus_event {
544         picture {
545             background_color = 255,255,255;

```

```

546         line_graphic grid;
547         x = 0; y = 0;
548
549         bitmap imagel;
550         x = 0; y = 0;
551     } Obj_test_pic;
552     target_button = 2,3;
553     port_code = 4;
554     code = "Remember Object Test";
555 };
556 } Obj_test_trial;
557
558 #Practice Remember Location test trial
559 trial {
560     trial_type = correct_response;
561     trial_duration = forever;
562     stimulus_event {
563         picture {
564             background_color = 255,255,255;
565             line_graphic grid;
566             x = 0; y = 0;
567
568             bitmap circle;
569             x = 0; y = 0;
570         } practice_Loc_test_pic;
571         target_button = 2,3;
572         port_code = 9;
573         code = "Practice Remember Location Test";
574     };
575 } practice_Loc_test_trial;
576
577 #Remember Location test trial
578 trial {
579     trial_type = correct_response;
580     trial_duration = forever;
581     stimulus_event {
582         picture {
583             background_color = 255,255,255;
584             line_graphic grid;
585             x = 0; y = 0;
586
587             bitmap circle;
588             x = 0; y = 0;
589         } Loc_test_pic;
590         target_button = 2,3;
591         port_code = 5;
592         code = "Remember Location Test";
593     };
594 } Loc_test_trial;
595
596 #Practice Remember Object and Location test trial
597 trial {
598     trial_type = correct_response;
599     trial_duration = forever;
600     stimulus_event {
601         picture {
602             background_color = 255,255,255;
603             line_graphic grid;
604             x = 0; y = 0;
605
606             bitmap imagel;
607             x = 0; y = 0;
608         } practice_Conjunction_test_pic;
609         target_button = 2,3;

```

```

610     port_code = 10;
611     code = "Practice Remember Obj+Loc Test";
612 };
613 } practice_Conjunction_test_trial;
614
615 #Remember Object and Location test trial
616 trial {
617     trial_type = correct_response;
618     trial_duration = forever;
619     stimulus_event {
620         picture {
621             background_color = 255,255,255;
622             line_graphic grid;
623             x = 0; y = 0;
624
625             bitmap image1;
626             x = 0; y = 0;
627         } Conjunction_test_pic;
628         target_button = 2,3;
629         port_code = 6;
630         code = "Remember Obj+Loc Test";
631     };
632 } Conjunction_test_trial;
633
634 #Inter-trial Blank Trial
635 trial {
636     trial_type = fixed;
637     trial_duration = 500;
638     picture {
639         box blank_box;
640         x = 0; y = 0;
641     };
642 } blank_trial;
643
644 #Rest between Part 1 and Part 2
645 trial {
646     trial_type = correct_response;
647     trial_duration = forever;
648     stimulus_event {
649         picture {
650             box blank_box;
651             x = 0; y = 0;
652
653             text {
654                 caption = "Rest";
655                 font_size = 44;
656             };
657             x = 0; y = 150;
658
659             text {
660                 caption = "Press the space bar to begin Part 2.";
661                 font_size = 36;
662             };
663             x = 0; y = -410;
664
665         } rest_pic;
666         target_button = 1;
667         code = "Rest";
668     } rest_event;
669 } rest_trial;
670
671 #Repeat Practice option for the experimenter
672 trial {
673     trial_type = correct_response;

```

```

674     trial_duration = forever;
675     stimulus_event {
676         picture {
677
678             text {
679                 caption = "Experimenter please press keys to continue.";
680                 font_size = 60;
681             };
682             x = 0; y = 110;
683
684         } repeat_practice_pic;
685         target_button = 4,5;
686         code = "Whether Repeat Practice";
687     };
688 } repeat_practice_trial;
689
690 #Practice Correct Feedback
691 trial {
692     trial_type = fixed;
693     trial_duration = 1500;
694     stimulus_event {
695         picture {
696
697             text {
698                 caption = "<b>CORRECT</b>";
699                 font_size = 60;
700                 font_color = 0,255,0;
701                 formatted_text = true;
702             };
703             x = 0; y = 110;
704
705         } feedback_correct_pic;
706         code = "Practice Correct Feedback";
707     };
708 } feedback_correct_trial;
709
710 #Practice Incorrect Feedback
711 trial {
712     trial_type = fixed;
713     trial_duration = 1500;
714     stimulus_event {
715         picture {
716
717             text {
718                 caption = "<b>INCORRECT</b>";
719                 font_size = 60;
720                 font_color = 255,0,0;
721                 formatted_text = true;
722             };
723             x = 0; y = 110;
724
725         } feedback_incorrect_pic;
726         code = "Practice Incorrect Feedback";
727     };
728 } feedback_incorrect_trial;
729
730 #####
731 #####
732 #House keeping parameters before beginning the experiment
733 #Begin PCL (Presentation Control Language)
734 begin_pcl;
735
736 #specifying output file

```

```

737 string logpath = logfile_directory;
738 string fn = logpath +logfile.subject()+"_ObjLoc.txt";
739 logfile.set_filename(logpath +logfile.subject() + "_ObjLoc_logfile.log");
740 output_file ofile1 = new output_file;
741 ofile1.open_append(fn);
742
743 #Setting column headers for output file
744 ofile1.print("Block\t"+"Trial\t"+"Cond\t"+"Studied\t"+"Resp\t"+"Corr\t"+"RT\n");
745
746 #Variables
747 int Block1 = 48; #This should be 48 after done testing
748 int Block2 = 24; #This should be 24 after done testing
749 int Studied_pic_count = 1;
750 int Non_studied_pic_count = 1;
751 int Loc;
752 int Resp;
753
754 response_data my_response;
755 response_data my_response_practice;
756 stimulus_data my_data;
757 stimulus_data my_data_practice;
758 int resp_button;
759 int correct_button;
760 int resp_button_practice;
761 int RT;
762 int corr_resp;
763 string condition;
764 int studied;
765
766 #Concatenate the two arrays of picture stimuli
767 pictures.append(pictures_b);
768 #Randomize the order of picture stimuli
769 pictures.shuffle();
770
771 #Instantiate empty bitmap arrays
772 array<bitmap>Studied_pics[240];
773 array<bitmap>Non_studied_pics[18]; # Reserve 18 pictures to be new stimuli that
will only be shown at testing
774
775 #Update array values with picture stimuli
776 loop int i = 1 until i > pictures.count() begin
777     if (i <= 240) then
778         Studied_pics[i] = pictures[i];
779     else
780         Non_studied_pics[i-240] = pictures[i];
781     end;
782     i = i + 1;
783 end;
784
785
786 #Trial Condition for Block1: 0 = Obj condition New, 1 = Obj condition Old, 2 =
Loc condition New, 3 = Loc condition Old
787 array<int>Block1_trial_type_index[48] = {0,0,0,0,0,0,0,0,0,0,0,0,0,0,1,1,1,1,1,1,1,1,1,
1,1,1,2,2,2,2,2,2,2,2,2,2,2,2,2,3,3,3,3,3,3,3,3,3,3,3,3};
788 Block1_trial_type_index.shuffle();
789
790 #Trial Condition for Block2: 0 = New Location, 1 = Old Location
791 #N.B.: conditions only concern location because the objects are always old; no
new object will be shown in Part 2.
792 array<int>Block2_trial_type_index[24] = {0,0,0,0,0,0,0,0,0,0,0,0,0,0,1,1,1,1,1,1,1,1,1,
1,1,1};
793 Block2_trial_type_index.shuffle();
794
795 #Trial Condition for practice1

```

```

796 array<int>Practice_trial_type_index[4] = {0,1,2,3};
797 array<int>Practice_correct_key[4] = {3,2,3,2};
798
799 #Trial Condition for practice2
800 array<int>Practice2_trial_type_index[2] = {0,1};
801 array<int>Practice2_correct_key[2] = {3,2};
802
803 #Possible object positions
804 array<double>object_position[8][2]={{-333.33,-333.33},{-333.33,0},{-333.33,333.33},
    },{0,-333.33},{0,333.33},{333.33,-333.33},{333.33,0},{333.33,333.33}};
805
806 #Object test position index
807 array<int>Obj_test_pos_index[3] = {1,2,3};
808
809 #Object test pic index
810 array<int>Obj_test_pic_index[3] = {1,2,3};
811
812 #####
813 #
814 # Beginning of experiment presentation
815 #
816 #####
817
818 #####
819 output_port port = output_port_manager.get_port( 1 );
820
821 #Instructions
822 #response_manager.set_port_output( false );
823 Instructions_trial.present();
824 Instructions_event.set_stimulus(instructions2);
825 Instructions_trial.present();
826 Instructions_event.set_stimulus(instructions3);
827 Instructions_trial.present();
828 Instructions_event.set_stimulus(instructions4);
829 Instructions_trial.present();
830 Instructions_event.set_stimulus(instructions5);
831 Instructions_trial.present();
832 Instructions_event.set_stimulus(instructions6);
833 Instructions_trial.present();
834 Instructions_event.set_stimulus(instructions7);
835 Instructions_trial.present();
836 Instructions_event.set_stimulus(instructions8);
837 Instructions_trial.present();
838 Instructions_event.set_stimulus(instructions9);
839 Instructions_trial.present();
840
841 #####
842 #Part 1 Practice - 4 trials: 1 Object-New, 1 Object-Old, 1 Location-New, 1
Location-Old
843 int Studied_pic_count_holder = Studied_pic_count;
844 int Non_studied_pic_count_holder = Non_studied_pic_count;
845
846 wait_interval (100);
847 port.send_code(22);
848 wait_interval (100);
849
850 loop int k = 0 until k > 0 begin
851     loop int j = 1 until j > 4 begin
852         blank_grid_trial.present();
853         object_position.shuffle();

```

```

854     loop int i = 1 until i > 3 begin
855         practice_pic.set_part(2, Studied_pics[Studied_pic_count]);
856         practice_pic.set_part_x(2, object_position[i][1]);
857         practice_pic.set_part_y(2, object_position[i][2]);
858         practice_study_trial.present();
859
860         wait_interval(20);
861
862         i = i + 1;
863         Studied_pic_count = Studied_pic_count + 1;
864     end;
865     wait_delay.present();
866
867     if (Practice_trial_type_index[j] == 0 || Practice_trial_type_index[j] == 1)
then #Obj condition
868         condition = "Obj";
869         Obj_Condition_Message.present();
870         if (Practice_trial_type_index[j] == 0) then #Not-studied Object trial
871             studied = 0;
872             correct_button = 3;
873             practice_Obj_test_pic.set_part(2,
Non_studied_pics[Non_studied_pic_count]);
874             practice_Obj_test_trial.present();
875             Non_studied_pic_count = Non_studied_pic_count + 1;
876         else #Studied Object trial
877             studied = 1;
878             correct_button = 2;
879             #Randomly choose one of the 3 pictures as the test stimulus
880             Obj_test_pic_index.shuffle();
881             practice_Obj_test_pic.set_part(2, Studied_pics[Studied_pic_count -
Obj_test_pic_index[1]]);
882             practice_Obj_test_trial.present();
883         end;
884     else #Loc condition
885         condition = "Loc";
886         Loc_Condition_Message.present();
887         if (Practice_trial_type_index[j] == 2) then #Not-studied Location trial
888             studied = 0;
889             correct_button = 3;
890             practice_Loc_test_pic.set_part_x(2, object_position[4][1]);
891             practice_Loc_test_pic.set_part_y(2, object_position[4][2]);
892             practice_Loc_test_trial.present();
893         else #Studied Location trial
894             studied = 1;
895             correct_button = 2;
896             Obj_test_pos_index.shuffle();
897             practice_Loc_test_pic.set_part_x(2,
object_position[Obj_test_pos_index[1]][1]);
898             practice_Loc_test_pic.set_part_y(2,
object_position[Obj_test_pos_index[1]][2]);
899             practice_Loc_test_trial.present();
900         end;
901     end;
902
903     wait_interval(20);
904
905     my_data = stimulus_manager.last_stimulus_data();
906     my_response = response_manager.last_response_data();
907     RT = my_data.reaction_time();
908     resp_button = my_data.button();
909
910     if (resp_button == Practice_correct_key[j]) then
911         feedback_correct_trial.present();
912     else

```

```

913         feedback_incorrect_trial.present();
914     end;
915
916     blank_trial.present();
917     j = j + 1;
918 end;
919
920 #Give experimenter option to repeat practice trials if necessary
921 repeat_practice_trial.present();
922
923 my_data = stimulus_manager.last_stimulus_data();
924 resp_button = my_data.button();
925
926 if (resp_button == 4) then
927     #Exit the outer for loop and proceed to Part 1
928     k = 1;
929 else
930     #Reset the counters so we re-use the practice trial pictures
931     Studied_pic_count = Studied_pic_count_holder;
932     Non_studied_pic_count = Non_studied_pic_count_holder;
933 end;
934
935 end;
936
937 wait_interval (100);
938 port.send_code(23);
939 wait_interval (100);
940 #Set the port output to be true - start writing to the output log
941 #response_manager.set_port_output( true );
942
943 #####
944 #####
945 #Part 1 Single Modality Trials
946 port.send_code(22);
947 wait_interval (100);
948
949 loop int j = 1 until j > Block1 begin
950     blank_grid_trial.present();
951     object_position.shuffle();
952     loop int i = 1 until i > 3 begin
953         pic.set_part(2, Studied_pics[Studied_pic_count]);
954         pic.set_part_x(2, object_position[i][1]);
955         pic.set_part_y(2, object_position[i][2]);
956         study_trial.present();
957
958         wait_interval(20);
959
960         i = i + 1;
961         Studied_pic_count = Studied_pic_count + 1;
962     end;
963
964     wait_delay.present();
965
966     if (Block1_trial_type_index[j] == 0 || Block1_trial_type_index[j] == 1) then
967 #Obj condition
968         condition = "Obj";
969         Obj_Condition_Message.present();
970         if (Block1_trial_type_index[j] == 0) then #Not-studied Object trial
971             studied = 0;
972             correct_button = 3;
973             Obj_test_pic.set_part(2, Non_studied_pics[Non_studied_pic_count]);
974             Obj_test_trial.present();
975             Non_studied_pic_count = Non_studied_pic_count + 1;
976         else #Studied Object trial

```



```

975         studied = 1;
976         correct_button = 2;
977         Obj_test_pic_index.shuffle();
978         Obj_test_pic.set_part(2, Studied_pics[Studied_pic_count -
Obj_test_pic_index[1]]);
979         Obj_test_trial.present();
980     end;
981     else #Loc condition
982         condition = "Loc";
983         Loc_Condition_Message.present();
984         if (Block1_trial_type_index[j] == 2) then #Not-studied Location trial
985             studied = 0;
986             correct_button = 3;
987             Loc_test_pic.set_part_x(2, object_position[4][1]);
988             Loc_test_pic.set_part_y(2, object_position[4][2]);
989             Loc_test_trial.present();
990         else #Studied Location trial
991             studied = 1;
992             correct_button = 2;
993             Obj_test_pos_index.shuffle();
994             Loc_test_pic.set_part_x(2, object_position[Obj_test_pos_index[1]][1]);
995             Loc_test_pic.set_part_y(2, object_position[Obj_test_pos_index[1]][2]);
996             Loc_test_trial.present();
997         end;
998     end;
999
1000     wait_interval(20);
1001
1002     #Saving data
1003     my_data = stimulus_manager.last_stimulus_data();
1004     my_response = response_manager.last_response_data();
1005     RT = my_data.reaction_time();
1006     # RT = my_response.reaction_time();??? - it should be my_data.reaciton_time(),
my_response doesn't have .reaction_time()
1007     resp_button = my_data.button();
1008
1009     if (resp_button == correct_button) then
1010         corr_resp = 1;
1011     else
1012         corr_resp = 0;
1013     end;
1014
1015     if (resp_button == 2) then
1016         Resp = 1;
1017     else
1018         Resp = 0;
1019     end;
1020
1021     ofile1.print(string(1) + "\t"); #Block number
1022     ofile1.print(string(j) + "\t"); #Trial number
1023     ofile1.print(condition + "\t"); #Condition type (Obj, Loc, or Both)
1024     ofile1.print(string(studied) + "\t"); #Was object studied (0=Not-studied;
1=Studied)
1025     ofile1.print(string(Resp) + "\t"); #Response (0=New; 1=Old)
1026     ofile1.print(string(corr_resp) + "\t"); # Was response correct (0=incorrect;
1=correct)
1027     ofile1.print(string(RT) + "\n"); #Response Time
1028
1029     blank_trial.present();
1030     j = j + 1;
1031 end;
1032
1033 wait_interval (100);
1034 port.send_code(23);

```

```

1035 wait_interval (100);
1036
1037 #Show the Rest trial, turn off port output during Rest
1038 #response_manager.set_port_output( false );
1039 rest_trial.present();
1040 #response_manager.set_port_output( true );
1041
1042 #####
1043 #Part 2 Instructions
1044
1045 #response_manager.set_port_output( false );
1046 Instructions_event.set_stimulus(instructions10);
1047 Instructions_trial.present();
1048 Instructions_event.set_stimulus(instructions11);
1049 Instructions_trial.present();
1050 Instructions_event.set_stimulus(instructions12);
1051 Instructions_trial.present();
1052 Instructions_event.set_stimulus(instructions13);
1053 Instructions_trial.present();
1054
1055 #####
1056 #Part 2 Practice - 2 trials: 1 New Location, 1 Old Location
1057 Studied_pic_count_holder = Studied_pic_count;
1058 Non_studied_pic_count_holder = Non_studied_pic_count;
1059
1060 wait_interval (100);
1061 port.send_code(22);
1062 wait_interval (100);
1063
1064 loop int k = 0 until k > 0 begin
1065     loop int j = 1 until j > 2 begin
1066         blank_grid_trial.present();
1067         object_position.shuffle();
1068         loop int i = 1 until i > 3 begin
1069             practice_pic.set_part(2, Studied_pics[Studied_pic_count]);
1070             practice_pic.set_part_x(2, object_position[i][1]);
1071             practice_pic.set_part_y(2, object_position[i][2]);
1072             practice_study_trial.present();
1073
1074             wait_interval(20);
1075
1076             i = i + 1;
1077             Studied_pic_count = Studied_pic_count + 1;
1078         end;
1079
1080         wait_delay.present();
1081         Conjunction_Condition_Message.present();
1082         condition = "Both";
1083
1084         if (Practice2_trial_type_index[j] == 0) then #Not-studied Location trial
1085             studied = 0;
1086             correct_button = 3;
1087             Obj_test_pic_index.shuffle();
1088             practice_Conjunction_test_pic.set_part(2, Studied_pics[Studied_pic_count
- Obj_test_pic_index[1]]);
1089             if (Obj_test_pic_index[1] == 3) then
1090                 Loc = random(2,3);
1091                 practice_Conjunction_test_pic.set_part_x(2, object_position[Loc][1]);
1092                 practice_Conjunction_test_pic.set_part_y(2, object_position[Loc][2]);
1093                 practice_Conjunction_test_trial.present();
1094             elseif (Obj_test_pic_index[1] == 2) then
1095                 int temp = random(1,2);

```

```

1096         if (temp == 1) then
1097             Loc = 1;
1098         else
1099             Loc = 3;
1100         end;
1101         practice_Conjunction_test_pic.set_part_x(2, object_position[Loc][1]);
1102         practice_Conjunction_test_pic.set_part_y(2, object_position[Loc][2]);
1103         practice_Conjunction_test_trial.present();
1104     else
1105         Loc = random(1,2);
1106         practice_Conjunction_test_pic.set_part_x(2, object_position[Loc][1]);
1107         practice_Conjunction_test_pic.set_part_y(2, object_position[Loc][2]);
1108         practice_Conjunction_test_trial.present();
1109     end
1110 else #Studied Location trial
1111     studied = 1;
1112     correct_button = 2;
1113     Obj_test_pic_index.shuffle();
1114     practice_Conjunction_test_pic.set_part(2, Studied_pics[Studied_pic_count
- Obj_test_pic_index[1]]);
1115
1116     if (Obj_test_pic_index[1] == 3) then
1117         practice_Conjunction_test_pic.set_part_x(2, object_position[1][1]);
1118         practice_Conjunction_test_pic.set_part_y(2, object_position[1][2]);
1119     elseif (Obj_test_pic_index[1] == 2) then
1120         practice_Conjunction_test_pic.set_part_x(2, object_position[2][1]);
1121         practice_Conjunction_test_pic.set_part_y(2, object_position[2][2]);
1122     else
1123         practice_Conjunction_test_pic.set_part_x(2, object_position[3][1]);
1124         practice_Conjunction_test_pic.set_part_y(2, object_position[3][2]);
1125     end;
1126     practice_Conjunction_test_trial.present();
1127 end;
1128
1129 wait_interval(20);
1130
1131 my_data = stimulus_manager.last_stimulus_data();
1132 my_response = response_manager.last_response_data();
1133 RT = my_data.reaction_time();
1134 resp_button = my_data.button();
1135
1136 if (resp_button == Practice2_correct_key[j]) then
1137     feedback_correct_trial.present();
1138 else
1139     feedback_incorrect_trial.present();
1140 end;
1141
1142 blank_trial.present();
1143 j = j + 1;
1144 end;
1145
1146 #Give experimenter option to repeat practice trials if necessary
1147 repeat_practice_trial.present();
1148
1149 my_data = stimulus_manager.last_stimulus_data();
1150 resp_button = my_data.button();
1151
1152 if (resp_button == 4) then
1153     #Exit the outer for loop and proceed to Part 1
1154     k = 1;
1155 else
1156     #Reset the counters so we re-use the practice trial pictures
1157     Studied_pic_count = Studied_pic_count_holder;
1158     Non_studied_pic_count = Non_studied_pic_count_holder;

```

```

1159     end;
1160
1161 end;
1162
1163 wait_interval (100);
1164 port.send_code(23);
1165 wait_interval (100);
1166 #Reset the port output to be true - continue writing to the output log
1167 #response_manager.set_port_output( true );
1168
1169 #####
1170 #####
1171 #Part 2 Conjunction Condition Trials
1172 port.send_code(22);
1173 wait_interval (100);
1174
1175 loop int j = 1; until j > Block2 begin
1176     blank_grid_trial.present();
1177     object_position.shuffle();
1178     loop int i = 1 until i > 3 begin
1179         pic.set_part(2, Studied_pics[Studied_pic_count]);
1180         pic.set_part_x(2, object_position[i][1]);
1181         pic.set_part_y(2, object_position[i][2]);
1182         study_trial.present();
1183
1184         wait_interval(20);
1185
1186         i = i + 1;
1187         Studied_pic_count = Studied_pic_count + 1;
1188     end;
1189
1190     wait_delay.present();
1191     Conjunction_Condition_Message.present();
1192     condition = "Both";
1193
1194     if (Block2_trial_type_index[j] == 0) then #Not-studied Location trial
1195         studied = 0;
1196         correct_button = 3;
1197         Obj_test_pic_index.shuffle();
1198         Conjunction_test_pic.set_part(2, Studied_pics[Studied_pic_count -
1199         Obj_test_pic_index[1]]);
1200         if (Obj_test_pic_index[1] == 3) then
1201             Loc = random(2,3);
1202             Conjunction_test_pic.set_part_x(2, object_position[Loc][1]);
1203             Conjunction_test_pic.set_part_y(2, object_position[Loc][2]);
1204             Conjunction_test_trial.present();
1205         elseif (Obj_test_pic_index[1] == 2) then
1206             int temp = random(1,2);
1207             if (temp == 1) then
1208                 Loc = 1;
1209             else
1210                 Loc = 3;
1211             end;
1212             Conjunction_test_pic.set_part_x(2, object_position[Loc][1]);
1213             Conjunction_test_pic.set_part_y(2, object_position[Loc][2]);
1214             Conjunction_test_trial.present();
1215         else
1216             Loc = random(1,2);
1217             Conjunction_test_pic.set_part_x(2, object_position[Loc][1]);
1218             Conjunction_test_pic.set_part_y(2, object_position[Loc][2]);
1219             Conjunction_test_trial.present();
1220         end
1221     else #Studied Location trial
1222         studied = 1;

```

```

1221     correct_button = 2;
1222     Obj_test_pic_index.shuffle();
1223     Conjunction_test_pic.set_part(2, Studied_pics[Studied_pic_count -
Obj_test_pic_index[1]]);
1224
1225     if (Obj_test_pic_index[1] == 3) then
1226         Conjunction_test_pic.set_part_x(2, object_position[1][1]);
1227         Conjunction_test_pic.set_part_y(2, object_position[1][2]);
1228     elseif (Obj_test_pic_index[1] == 2) then
1229         Conjunction_test_pic.set_part_x(2, object_position[2][1]);
1230         Conjunction_test_pic.set_part_y(2, object_position[2][2]);
1231     else
1232         Conjunction_test_pic.set_part_x(2, object_position[3][1]);
1233         Conjunction_test_pic.set_part_y(2, object_position[3][2]);
1234     end;
1235     Conjunction_test_trial.present();
1236 end;
1237
1238 wait_interval(20);
1239
1240 #Saving data
1241 my_data = stimulus_manager.last_stimulus_data();
1242 my_response = response_manager.last_response_data();
1243 RT = my_data.reaction_time();
1244 resp_button = my_data.button();
1245
1246 if (resp_button == correct_button) then
1247     corr_resp = 1;
1248 else
1249     corr_resp = 0;
1250 end;
1251
1252 if (resp_button == 2) then
1253     Resp = 1;
1254 else
1255     Resp = 0;
1256 end;
1257
1258 ofile1.print(string(2) + "\t"); #Block number
1259 ofile1.print(string(j) + "\t"); #Trial number
1260 ofile1.print(condition + "\t"); #Condition type (Obj, Loc, or Both)
1261 ofile1.print(string(studied) + "\t"); #Was object studied (0=Not-studied;
1=Studied)
1262 ofile1.print(string(Resp) + "\t"); #Response (0=New; 1=Old)
1263 ofile1.print(string(corr_resp) + "\t"); # Was response correct (0=incorrect;
1=correct)
1264 ofile1.print(string(RT) + "\n"); #Response Time
1265
1266 blank_trial.present();
1267 j = j + 1;
1268 end;
1269
1270 #####
1271 #End of experiment presentation
1272 wait_interval (100);
1273 port.send_code(23);
1274 wait_interval (100);
1275
1276 #close output file
1277 ofile1.close();

```